shouldfix open (sfb)

REF. 8410

Arm-Schulter
Fixationsbandage

Arm-shoulder immobiliser

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ENGLISH

Please read the instructions carefully and thoroughly

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/EEC and 2007/47/EC. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro® Industries B.V.

MAINTENANCE

- Do not bleach No chemical cleaning
- Do not iron Do not tumble-dry
- Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. The device should not be worn in the vicinity of open flames.

WARNINGS

The product may only be used in the event of medical indication or when prescribed by a physician and must be applied by an orthopaedic technician, the competent authority figure for the application as well as for information on the safe use, according to individual needs. To ensure the effectiveness, tolerability, and correct operation, it is necessary that the device is put on with the utmost care. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. Use of the product is recommended for a single patient. With hypersensitive patients, direct contact with skin may cause redness or irritation. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician.

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Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland

Tel. 0041 91 8220088 · Fax 0041 91 8220089

info@orthoservice.com · www.orthoservice.com

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Tel. 0049 (0)7221 9719780 · Fax 0049 (0)7221 9719782

info@orthoservice.de · www.orthoservice.de

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ENGLISH

CHARACTERISTICS AND MATERIALS

- Structure in soft bonded fabric with the outer layer in brushed nylon to facilitate the Velcro® fastening
- Shoulder strap with side-release buckle closure

INDICATIONS

- Simple humeral fractures
- Glenohumeral luxation
- Distorsions of the shoulder
- Postoperative care of humeral fractures and interventions due to dislocations of the shoulder
- Postoperative care for tendon sutures of the long head of the biceps muscle

CONTRAINDICATIONS

Currently no known

APPLICATION

- 1 Pass the waist belt (A) around the waist and use the Velcro® closure to tension firmly and close.
- 2 Bring the arm cuff (B) that is sewn to the waist belt up to the height of the affected arm and then secure by closing the appropriate Velcro® closure.
- 3 After inserting the arm, secure the lower arm support brace (C) to the waist belt (A).
- 4 On the rear, secure the closure strap that encircles the shoulder (D) to the waist belt (A) using the Velcro® closure. Pass the strap (D) over the shoulder opposite of the affected one and secure to the lower arm support brace with the clip buckle closure (E).
- 5 Tighten the strap (D) appropriately. In doing so, ensure that the protective pad for the neck (F) is positioned correctly.
- 6 Fix the upper part of the lower arm support brace (C) to the waist belt (A) using the Velcro® securing strap (G).

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